


Acquire raw projections
(sinogram)


napari



Load stack

1

napari-tomodl



Select image layer

2

Half-rotation scans

3

Manual or automatic
axis alignment

4

Data pre-processing
image size, flat-field correction

5

Reconstruction methods

6

FBP TwIST U-Net ToMoDL

Reconstruction settings

filter, smoothing level, clip to circle

7

Reconstruct full volume
or specific slice(s) in # batches

8

Choose rotation axis

9

Volume post-processing
invert colours, 16-bit conversion

10